A Parametric Study of Electron Backscatter Diffraction based Grain Size Measurements

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Angenommen: 13. November 2009

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Abstract

Polycrystalline microstructures can be well imaged using maps reconstructed from orientation data collected by electron backscatter diffraction (EBSD) in the scanning electron microscope (SEM). These maps are very helpful for delineating grain boundaries in such microstructures and are thus well suited for grain size analysis. However, for these analyses to be as accurate as possible it is important to understand the effect of various parameters relevant to grain size measurements via EBSD in order to achieve accurate estimations of grain size. In particular, this work addresses the sensitivity of the grain tolerance angle and minimum number of measurement points per grain to the end results as well as considers other relevant parameters. Guidelines are proposed for making appropriate value selections for these parameters.

Introduction

Electron backscatter diffraction (EBSD) in the scanning electron microscope (SEM) is a powerful tool for characterizing the crystallographic orientation aspects of polycrystalline materials. Although not always recognized as such, grain size is an orientation related microstructural feature.; as individual grains have unique crystallographic orientation relative to one-another. Automated EBSD scans collect orientation data from discrete points on a regular grid (often

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termed orientation imaging microscopy or OIM [1]). At each point on the grid, an EBSD pattern is collected and the orientation automatically determined from the pattern through a process called indexing. The orientation (along with other parameters such as a measure of indexing confidence) is recorded. After the orientation data is collected, neighboring points on the grid with similar orientation can be grouped together into "grains" [2, 3]. These grain groups typically correspond to the grains observed in conventional metallography as shown in Fig. 1. Early EBSD systems were relatively slow and would require overnight scans to collect enough scan data to reconstruct the polycrystalline microstructure from the scan data points. However, the speed of modern EBSD systems has increased dramatically in recent years. Thus, the characterization of grain size in polycrystals using EBSD is now much more practical. EBSD is particularly helpful in characterizing aluminum alloys where traditional metallographic analysis can be difficult due to the lack of contrast in optical micrographs.

Despite the advances in EBSD data acquisition, analysis of EBSD data for grain size measurement remains somewhat difficult due to a lack of standard definition for several key parameters [4-9]. The most critical parameters involved in constructing grains from the EBSD data are the resolution of the scan and the tolerance angle required to define orientation similarity. The selection of the tolerance angle defining the orientation similarity has a strong effect on the size of the grains resulting from the grain grouping algorithm. The choice of grain tolerance angle is particularly problematic in deformed materials [10, 11]. In this study, we limit the analysis to recrystallized materials.

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The choice of values for several key parameters is based on comparisons with ASTM standard E1382 – Standard Test Methods for Determining Average Grain Size Using Semiautomatic/Automatic Image Analysis. The parameters considered are (1) the tolerance angle, (2) the step size between points in the grid, (3) the minimum number of measurements (or pixels) for a group of similarly oriented grid points to be considered a “grain” and (4) the overall size of the scan grid. This work shows the sensitivity of the grain size measurements to these parameters and proposes guidelines for selecting appropriate values. In addition errors unique to EBSD measurements will be discussed.

Data

Two different sets of data are used in this work to explore the effect of the various parameters on the resulting grain size measurements. One data set shows a relatively equiaxed grain size structure from an as-deposited aluminum thin film and the other data set exhibits an elongated structure from a polished section of a recrystallized bulk 6000 series aluminum alloy sample. The data set obtained from the equiaxed structure contains just over four million data points. The data set from the elongated structure contains nearly 580 thousand data points. The data points in the scan grids are oft times referred to as pixels particularly when considering maps constructed from the data sets.

Maps reconstructed from these two datasets are shown in Fig. 1. The colors in these maps are essentially random but chosen so that no two neighboring grains are shaded using the same color so that the grains can be delineated by color. The color maps are overlaid on gray scale maps based on image quality (IQ). The IQ is a measure of the quality of the diffraction pattern at each point in the OIM scan. This is mapped to a gray scale to produce the map shown. The

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
important thing to note in the IQ map, is that grain boundaries are delineated by intensity as opposed to orientation as is done in the color overlay. When the beam is positioned on a near a grain boundary in the microstructure, the diffraction volume will include the crystal lattices from both sides of the boundary. This results in an EBSD pattern that is a composite of the individual patterns from the two grains on either side of the boundary. These patterns typically produce much lower IQ values. It should be noted that no post processing cleanup was applied to the data for the maps shown in Fig. 1.

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


Figs. 1a and b. Color grain maps overlaid on image quality gray scale maps from OIM data on (a) a recrystallized bulk 6000 series aluminum alloy sample and (b) an aluminum thin film

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Grain size analysis of the equiaxed data set was performed using an equivalent area approach. In this approach, the number of points in each grain is analyzed. The area encompassed by a grain is proportional to the number of data points contained within the grain. For a data set collected on a square the grid the area of the grain would simply be then number of data points multiplied by the square of the step size between scan grid points. For the hexagonal scan grids used in these two data sets, the area of a grain is the number of points in the grain multiplied by the square of the

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Figs. 2a to c. Color grain maps overlaid on image quality maps for varying tolerance angles of (a) 1°, (b) 5° and (c) 15°. The sample is an as-deposited aluminum thin film.

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

step size and in turn multiplied by $\sqrt{3}/2$. A circle of equivalent area is constructed and an equivalent diameter can then be determined.

A linear intercept method was used to analyze the grain size in the elongated structure. Both vertical and linear intercepts were used in order to characterize the elongated structure.

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Figs. 3a and b. Average grain size plotted as a function of the grain tolerance angle for (a) the elongated structure and (b) the equiaxed structure

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Grain Tolerance Angle


To use EBSD data for grain size analysis, a tolerance angle defining whether two neighboring points belong to the same grain must be selected. The grain size results can be very sensitive to the choice of grain tolerance angle. This is evident in Fig. 2, where different tolerance angles are used in the grain reconstruction algorithm. These results are shown overlaid on an image quality (IQ) map.

It is evident from Fig. 2 that the “grain size” increases as the tolerance angle increases. Essentially, grains separated by low angle boundaries coalesce into single grains at larger tolerance angles. Even at 1° tolerance angle there appear to be grain boundaries in the IQ map that are not captured in the grain grouping algorithm. The choice of grain tolerance angle is; thus, really a decision of whether such features represent one grain or two [12]. The decision as to the correct choice of tolerance angle can be difficult. A greater understanding of the role of these small angle boundaries in property models where grain size is assumed to play a role would aid in the decision making process. In

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Figs. 4a to d. Grain maps reconstructed from EBSD scan data from recrystallized aluminum at step sizes of (a) 1, (b) 2, (c) 4 and (d) 8 μm . The colors in these maps are randomly selected to delineate individual grains. White points are grains with too few pixels to be considered as grains. Black points are non-indexed data points.

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
Fig. 5. Average grain size plotted as a function of step size for the elongated structure.

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order to better quantify the effect of the choice of tolerance angle Fig. 3 plots the measured average grain size as a function of tolerance angle. For the elongated data set, the average grain size is relatively insensitive to the choice of grain tolerance angles from 1° to 15° whereas the thin film data is more sensitive to the choice of grain tolerance angle. In this case, the curve appears to plateau at values between 4° to 8°. For both data sets, a choice of 5° for the grain tolerance angle appears appropriate. This is fortuitous, as this is a default value used in commercial EBSD packages which is seldom changed by users.

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Step Size

In order to explore the effects of scan step size on the grain size analysis, the scan data for the elongated was coarsened by removing every other row and every other column of data. This was repeated to create four scans with 1, 2, 4, and 8 µm step sizes covering the same scan area as shown in Fig. 1. Fig. 4 shows close-up views of the bottom left corner of the scan for all four step sizes.

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Fig. 5 shows the effect of increasing step size on the measured average grain size. The results

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
Fig. 6. Average grain size plotted as a function of the minimum grain size for the elongated structure

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suggest that the scan density associated with a step size of 4 µm is an acceptable minimum. This results in a scan density of approximately 71 points per average grain.

Minimum Grain Size

The step size is closely intertwined with the choice of minimum grain size. The minimum grain size is more precisely the minimum number of pixels before a grain group is considered a valid grain.

From the ASTM standard, the minimum sized grain that should be considered would possess a grain diameter of approximately 11 pixels. (This was derived by assuming a measurement field on a 512 x 512 pixel array on a 225 x 225 m² monitor. The stated minimum grain diameter is 5 mm, thus (5 mm x 512 pixels)/225 mm = 11.4 pixels.) If a circular grain shape is assumed, then the equivalent area would be 102 pixels. This is substantially larger than the values of 2-5 pixels typically used for EBSD grain size measurements. The small values are used since, two or three neighboring points were indexed with the same (or very similar) orientations during an automated scan then those points are assumed to be correctly indexed, particularly if they fall on multiple rows. To find a meaningful value the average grain size was determined using increasing minimum number of points as shown in Fig. 6 for the scans on the elongated material with 1 and 2 µm step sizes.

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
Fig. 7. Average grain size plotted as a function of step size for the equiaxed structure for a fixed minimum grain size of 100 pixels and for a minimum grain size that scales with step size (100, 50, 25, 13 nm).

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For both the 1 and 2 µm step size scans, the average grain size is relatively constant for minimum grain sizes between 3 and 20 pixels. Beyond 20 pixels, the scan resolution begins to play a large role in the measured grain areas. As would be expected, the effect of minimum grain size exhibits itself at smaller values for the coarser scan than the finer scan. For smaller minimum grain sizes the non-indexed or mis-indexed points will have a greater impact on results. The impact will be greater in data sets with a high-fraction of suspect points. The effects of data quality will be discussed in more detail in a later section.

For the ASTM standard, the ratio of the smallest grain size to the average grain size is approximately 100 to 5000. A similar ratio can be achieved if a minimum grain size of 10 pixels is used and the average pixels per grain is 500. This concept was used in creating Fig. 6. In this plot, the minimum grain size was scaled with increasing step size. This approach was also followed in Fig. 7 for the data from the equiaxed structure. In this plot, the error bars are generated by assuming a value of $\pm 10\%$ difference in the minimum grain size. The red curve was created without any scaling assuming a minimum grain size of 100 pixels regardless of the scan step size, whereas the blue curve was generated with a minimum grain size of 100 pixels for the full 4 million point scan and then scaling by a factor

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of four for a coarsening of the dataset by two in both the horizontal and vertical directions.

For the most part, maintaining an appropriate ratio between the minimum grain size and the average number of pixels per grain is more critical to achieving appropriate grain size results than actual absolute values.

Measurement Area

ASTM E1382 recommends that the scan area contains approximately 50 grains. A typical 512×512 pixel image contains 262,144 pixels. This is a considerable number of measurements for an EBSD scan. Although modern EBSD systems are capable of scan speeds of over 200 measurements per second with reliable indexing on field emission gun SEMs, speeds of 50 points per second are more typical for the installed base of systems. If 5 fields are to be measured then this represents a measurement time of 1.8 h with modern EBSD systems and 7.3 h for older systems. Grain boundaries are quantitatively defined in automated EBSD scans, as opposed to simple image intensity deviations in optical measurements. Thus, image processing techniques such as skeletonization or watershed segmentation are not typically used in grain size measurements performed using EBSD. For these reasons, a sparser measurement array is justified for accurate measurements of grain size using EBSD than used in optical measurements using digitized images. Note also that a 512×512 scan area enclosing 50 grains means that the average grain would contain over 5000 pixels.

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Other Sources of Error


Data Quality

In any EBSD scan, a certain fraction of the measurement points are either zero-solutions as the corresponding pattern could not be indexed or the orientation measured is incorrect. Various parameters can be used to judge the quality of the automatically indexed orientation data. One metric is the confidence index (CI) [13, 14]. Fig. 8 shows a map obtained from the equiaxed data with data with CI values less than 0.2 excluded along with a plot of the average grain size as a function of the CI filtering value. The CI filter preferentially filters out points near grain boundaries. As previously noted, points near grain boundaries typically produce patterns which are a superposition of the patterns from the two grains on either side of the boundary even though they may be indexed correctly. Thus, filtering the data by CI will tend to decrease the average grain size as points near grain boundaries will be excluded from the analysis even when they are indexed correctly.

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
Figs. 8a and b. (a) A grain map for the equiaxed structure with data points with CI less than 0.2 excluded from the map accompanied by (b) a plot of the average grain size versus the magnitude of the CI filter

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Figs. 9a to k. Grain maps obtained by re-scanning a set of saved EBSD patterns on steel where a percentage of the pixels in the patterns have been assigned random intensity – (a) 0%, (b) 20%, (c) 40%, (d) 60% and (e) 80%; (f-j) corresponding grain maps after applying a cleanup algorithm; (k) an accompanying plot of the grain size distribution of the raw and cleaned up data

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Various clean-up routines have been devised to differentiate points with low CI but are most likely indexed correctly from mis-indexed and non-indexed points and then to “clean-up” these points. This is typically done by re-assigning the orientation of the non-indexed or suspect points with the orientation of one of the neighboring points on the scan grid [15, 16]. Wright et al. [16] have shown the effect of clean up routines on the measured texture for a set of data from steel. The patterns at each point in the dataset were saved. The patterns were then re-indexed after randomly changing a fraction of the points in the pattern to simulate adding noise to the patterns. Thus, datasets were created from the noisy patterns with different amounts of added noise as shown in Fig. 9. The data were then cleaned and the grain size tracked and compared to the data obtained without any added noise (labeled 0% noise in Fig. 9). In all cases, the cleaned data led to distribution curves more closely aligned with the distributions obtained with the 0% noise data. However, other instances have been observed when the cleanup routines lead to faulty results. Careful cleanup procedures generally improve the grain size results. Of course, any cleanup artifacts can be avoided by starting with the best data possible.

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
Fig. 10. The error in estimating the size of the circle by a regular grid of pixels plotted as a function of the diameter of the circle in grid pixels accompanied by corresponding power law trend lines and square and hexagonal grid examples for a 8 pixel wide circle

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Equivalent Circle Method

The ability to accurately approximate a circle with a grid is a function of the resolution of the grid. Thus, for smaller grains, the equivalent circle approach to grain size measurements becomes less accurate. Fig. 10 shows the errors associated with the approximation of the area of a circle as a function of the number of points in a square and a hexagonal grid. However, if very many grains are measured the errors arising from this effect should be small, particularly if the average grain size analysis follows the proposed guidelines set forth below.

Tilt Effects

As Nolze [17] has pointed out, the high specimen tilts required for EBSD (typically 70°) can lead to various distortions. In fact, the accuracy of the measured tilt can be an issue. It is difficult during sample preparation to keep both sides of a polished sample exactly in parallel which can lead to inaccuracies depending on how the

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
Fig. 11. A plot showing the effects of tilt inaccuracies on data obtained from the elongated structures

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sample is mounted in the microscope. Fig. 11 shows the impact of the tilt on the grain size results obtained from the elongated structure. As would be expected, the impact in the horizontal direction is minor. However, the tilt in the vertical direction clearly affects the grain size analysis. One way to make sure the sample is at assumed position is to rotate the sample 90° about the sample normal and to rescan the same area. If the results are identical then the tilt of the sample is as assumed. This is reflected in Fig. 11 for the two clear larger markers at 71°.

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Conclusions

While the conventional standard is to measure 50 grains per field, a larger number of approximately 300 is recommended to improve the statistical sampling and reduce edge effects. For the parameters described here the following values are recommended for practical grain size analysis via EBSD on recrystallized materials:

- 150,000 scan points
- Minimum grain size of 10 pixels,
- 500 pixels for the average grain (150,000/300),
- Tolerance angle of 5°.

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This study has been performed using hexagonal scan grids; however, experience has shown that square scan measurement grids lead to similar conclusions. While these guidelines can help lead to reproducible grain size results, it is important to also perform scans so as to produce data of the best quality possible with high indexing success rates. Proper sample preparation and indexing procedures should be followed to avoid introducing any errors or biases into the data and care exercised in any post processing to clean-up the data.

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Acknowledgements

The author gratefully acknowledges helpful discussions with Scott Sitzman of Oxford-HKL, Matt Nowell of EDAX-TSL, David Field of Washington State University and Peter Quedstedt of the National Physical Laboratory in the UK. In addition, the author acknowledges Hasso Weiland and Ed Llewellyn of Alcoa for supplying the aluminum samples and Matt Nowell of EDAX-TSL for performing the OIM scans.

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stuart.wright@ametec.com

Stuart I. Wright

developed a fully automated EBSD system for his PhD at Yale University. Postdoc. at Los Alamos National Lab in texture characterization. Founding member of TSL (now part of EDAX) to commercialize the automated EBSD technique. Continues both technique and applications development on EBSD at EDAX in Salt Lake City, Utah.

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